

# The Dermillume Pro Series!

**Red and Blue** for Acne Treatments,  
**Red + NIR** for Anti-Aging Services

**Offer both Acne Treatments and Anti-aging Treatments  
with One Phototherapy Device**

- Use in **Red + Blue** mode for a drug-free, effective Acne treatment
- Use in **Blue-only** mode to activate generic ALA (aminolevulinic acid) for medical Acne treatments
- Use in **Red + Near Infrared** mode for Anti-aging treatments
- Use in **Red + Near Infrared** mode to enhance collagen treatments

**Visible results in one week**



**Acne  
& Skin  
Healing**



The Pro1500 Covers 6" x 13"  
20' treatment Red, NIR + Blue.



The Pro APL covers 2.5" x 3"  
90 sec treatment Red +Blue.

**Dermillume**  
Precision Light Therapy

**Safe - Fast - Effective**

## The Dermillume Advantage

- FDA approved for market for Acne treatments
- Field-tested anti-aging treatments
- 3 light programs
- Each device certified for power delivery to skin
- Does not touch skin
- Easy to use – minimal training

**Care Electronics, Inc.**  
303-444-2273 or 888-444-8284

Detailed specifications on reverse or at:

[www.dermillume.com](http://www.dermillume.com)

## Specifications for the Dermillume Pro Series

Improve Your Acne and Anti-aging Treatments by adding Proven Dermillume Pro Phototherapy to your existing services

- ☞ **Blue Light** kills acne bacteria deep in the skin where topical treatments cannot reach
- ☞ **Red + Near Infrared light smooths** wrinkles, tightens pores and smooths weathered or damaged skin by increasing the number of new skin cells in the damaged area to normal levels
- ☞ **Blue and Red Light** provide superior drug-free Acne treatments and Lesion Healing

### Dermillume Pro 1500 is Field Proven and Salon Tested

#### Power Delivered at the Skin

Pure 414 nm Blue light at 31 Joules/cm<sup>2</sup>  
Pure 660 nm Red at 85 Joules/cm<sup>2</sup>  
Pure 940 nm NIR at 10 Joules/cm<sup>2</sup>  
126 total Joules/cm<sup>2</sup>

#### Three Light Programs

Blue-only, 8', for ALA activation  
Blue/Red, 20', for drug-free acne healing  
Red + NIR, 20', for anti-aging treatments

Hands Free operation at 4 inches from skin  
20-minute treatment time – automatically controlled

Powered by plug In DC supply 18 Volts @ 3.5 Amps, using 120 volt wall socket

Warranty: 2-Year Parts and Labor

---

### Dermillume Pro Handheld APL is Field Proven and Salon Tested

#### Power Delivered at the Skin

Pure 414 nm Blue light at 18 Joules/cm<sup>2</sup>  
Pure 660 nm Red at 42 Joules/cm<sup>2</sup>  
60 total Joules/cm<sup>2</sup>

#### Three Light Programs

Blue-only, 45 sec, for ALA activation  
Blue/Red, 90 sec, for drug-free acne healing  
Red-only, 90 sec, for anti-aging treatments

Note: The Handheld does not include Near Infrared.

Warranty: 1-Year Parts and Labor