

The Dermillume Pro Hand Held APL!

Red and Blue for Acne Treatments, Red for Rejuvenation Services

Improve Acne Treatments with affordable phototherapy

- ☞ **Blue Light** kills bacteria deep in the skin where topical treatments can't reach
- ☞ **Red light** promotes healing of lesions by increasing the number of new skin cells in the damaged area

Activates generic ALA
(aminolevulonic acid) in 45 sec.



Dermillume Pro HHAPL in the Blue-only program

Offer both Acne Treatments and Anti-aging Treatments with One Device

- ☞ Use in **Red + Blue** mode for a drug-free, effective Acne treatment
 - ☞ Use in **Blue-only** mode with generic ALA (aminolevulonic acid) in dermatological settings
 - ☞ Use in **Red-only** mode for rejuvenation treatments
 - ☞ Use in **Red-only** mode to enhance collagen treatments – visible results in one week
-



Dermillume Pro HHIPL in the Red-only program



The Dermillume Advantage

- ☞ **60 Joules/cm² at skin in 90 seconds**
- ☞ **Treats 2.5" x 3" area**
- ☞ **Three programs deliver Blue-only, Red only or Red + Blue modes**
- ☞ **Use to activate Levulan™**
- ☞ **Does not touch skin**
- ☞ **Easy to use – minimal training**
- ☞ **Cordless Design**

In stock at Care Electronics, Inc.
(303) 444-2273 or (888) 444-8284

Optical Specifications for the Pro Hand Held APL™:

Blue wavelength – 414 nm

Red wavelength – 660 nm

- ☞ **Delivers 9.0 Joules/ cm² Blue-only at the skin in 45 seconds for ALA activation**
- ☞ **Delivers 42.0 Joules/ cm² Red-only at the skin in 90 seconds for Rejuvenation treatments**
- ☞ **Delivers 60.0 Joules/ cm² Red + Blue at the skin in 90 seconds for drug-free acne treatments**

More information at www.dermillume.com

Levulan is a trademarked product
of Dusa Pharmaceuticals