

The Dermillume Pro Hand Held IPL!

Red and Blue for Acne Treatments, Red for Rejuvenation Services

Improve Acne Treatments with affordable phototherapy

- ☞ **Blue Light** kills bacteria deep in the skin where topical treatments can't reach
- ☞ **Red light** promotes healing of lesions by increasing the number of new skin cells in the damaged area

Activates generic ALA
(aminolevulonic acid) in 45 sec.



Dermillume Pro HH IPL in the Blue-only program

Offer both Acne Treatments and Anti-aging Treatments with One Device

- ☞ Use in **Red + Blue** mode for a *drug-free, effective Acne treatment*
 - ☞ Use in **Blue-only** mode with *generic ALA (aminolevulonic acid)* in dermatological settings
 - ☞ Use in **Red-only** mode for *rejuvenation treatments*
 - ☞ Use in **Red-only** mode to *enhance collagen treatments* – visible results in one week
-



Dermillume Pro HH IPL in the Red-only program



The Dermillume Advantage

- ☞ **60 Joules/cm² at skin in 90 seconds**
- ☞ **Treats 2.5" x 3" area**
- ☞ **Three programs deliver Blue-only, Red only or Red + Blue modes**
- ☞ **Use to activate Levulan™**
- ☞ **Does not touch skin**
- ☞ **Easy to use – minimal training**
- ☞ **Cordless Design**

Care Electronics, Inc
303-444-2273 or 888-444-8287

Optical Specifications for the Pro Hand Held IPL™:

Blue wavelength – 414 nm

Red wavelength – 660 nm

- ☞ **Delivers 9.0 Joules/ cm² Blue-only at the skin in 45 seconds for ALA activation**
- ☞ **Delivers 42.0 Joules/ cm² Red-only at the skin in 90 seconds for Rejuvenation treatments**
- ☞ **Delivers 60.0 Joules/ cm² Red + Blue at the skin in 90 seconds for drug-free acne treatments**

More information at www.dermillume.com

Levulan is a trademarked product
of Dusa Pharmaceuticals