

# The Dermillume Salon Healing Red!

With Red and Near Infrared for a Superior Anti-aging Treatment

Improve Your Anti-Aging Treatments  
with affordable phototherapy

*So Much Better than a Wand!*

- ☞ **Red light** increases the number of important skin cells in damaged areas  
Fills wrinkles and fine lines naturally  
Reduces pore size
- ☞ **Red light** increases collagen in the skin  
For smoother, more elastic skin
- ☞ **Red light** increases the replacement rate of skin cells  
Naturally reverse sun damage without burning  
Reduce visibility of scar and stretch marks
- ☞ **Near Infrared light (NIR)** improves blood circulation under the skin for better tissue oxygenation  
Improves healing time of skin damage  
Enhances efficiency of red light



Accelerate lightening and wrinkle reduction treatments with Dermillume Healing Red

- ☞ Visible improvement in one week when used with collagen peptide treatments
- ☞ 4-6 weekly treatments give best results, with a monthly maintenance schedule

**Important!** - Only devices approved for market by the FDA are acceptable to your insurer.

In stock at Care Electronics, Inc.  
(303) 444-2273 or (888) 444-8284



- ☞ Two colors in one 20' treatment
- ☞ Delivers 105 Joules/cm<sup>2</sup> Red + NIR light in 20' place 4 " from skin
- ☞ Uniform, hands-free treatment
- ☞ 11" x 7" – full face or neck
- ☞ Automatically turns off
- ☞ Red-only modality for clients with Rosacea
- ☞ Easy to use
- ☞ Suitable for clients seated in chairs or reclining on treatment couches
- ☞ Compact storage

More information at [www.dermillume.com](http://www.dermillume.com)

Optical Specifications for the Salon Healing Red:

Red Wavelength - 660 nm  
Red Spectral Power Output - 6321 mWatt/cm<sup>2</sup>  
Delivers 62 Joules/cm<sup>2</sup> at the skin

NIR Wavelength - 940 nm  
NIR Spectral Power Output - 3272 mWatt/cm<sup>2</sup>  
Delivers 43 Joules/cm<sup>2</sup> at the skin